

QUICK STEPS TO A HEALTHIER YOU!

Step 1 Setting up

Set the batteries in your CHIEF weighing scale's Battery Slot

Download 'CHIEF UMCH' app to your mobile device

Step 2 Sign up & Link

Create your CHIEF Login ID and user profile




Link your CHIEF Body Scale to your mobile device

Step 3 Weigh & Sync

Congratulations!

Easier Weight Management for a Healthier YOU!

CHIEF Smart Body Fat User Guide

<p>Step 1 Setting up</p> <ul style="list-style-type: none"> • Set the required batteries in your CHIEF weighing scale (location on the back of the weighing scale) • For quick access, scan the QR Code to download the app • OR for Android users, download 'CHIEF UMCH' app from Google Play Store to your mobile device • For iOS users, download 'iWellness' app from App Store to your mobile device 	
<p>Step 2 Sign up & Link</p> <p>Sign up for your FREE CHIEF Login ID, create your profile and Link your CHIEF weighing scale to your mobile device. For existing users, Login using your email (CHIEF ID)</p> <ul style="list-style-type: none"> • Select 'Menu', followed by 'Link Device' at CHIEF app • Ensure Bluetooth on your mobile device is turned ON • Step on the weighing scale with bare feet to activate the weighing scale • Select 'Start Pairing' to Link CHIEF weighing scale to your mobile device 	
<p>Step 3 Weigh & Sync</p> <ul style="list-style-type: none"> • Step on your CHIEF weighing scale with bare feet • Keep standing on the weighing scale until your weight reading is displayed on your CHIEF app • Press 'Save' to save data to CHIEF app • Easier weight management for a Healthier YOU starts NOW! 	
<p>Follow us on ...</p> <ul style="list-style-type: none"> • Facebook: CHIEF UMCH • Youtube: CHIEF UMCH • WeChat: CHIEF_UMCH • WhatsApp: +6 018 218 1328 • CHIEF Smart Activity & Sleep Tracker Video Guide https://youtu.be/XxWbikmczaY 	